

For Research Use Only

# SLC29A3 Recombinant antibody, PBS Only

Catalog Number:83036-3-PBS



## Basic Information

Catalog Number:

83036-3-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18563

GenBank Accession Number:

BC120996

GeneID (NCBI):

55315

UNIPROT ID:

Q9BZD2

Full Name:

solute carrier family 29 (nucleoside transporters), member 3

Calculated MW:

475 aa, 52 kDa

Observed MW:

51 kDa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

240123A8

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)


W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83036-3-RR (SLC29A3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83036-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83036-3-RR against Human SLC29A3 were performed. The affinity constant is 0.824 nM.