

For Research Use Only

MOG Recombinant antibody, PBS Only (Detector)

Catalog Number: 83063-7-PBS



Basic Information

Catalog Number: 83063-7-PBS	GenBank Accession Number: N/A	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4340	CloneNo.: 250235F12
Source: Rabbit	UNIPROT ID: Q16653-1	
Isotype: IgG	Full Name: myelin oligodendrocyte glycoprotein	
Immunogen Catalog Number: EG2878	Calculated MW: 24 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

83063-7-PBS targets MOG as part of a matched antibody pair:

MP00110-4: 83063-8-PBS capture and 83063-7-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.



Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00110-4, Human MOG Recombinant Matched Antibody Pair - PBS only. 83063-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2878. 83063-7-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.

The mean MOG concentration was determined to be 158.5 ng/mL in mouse brain extract based on a 6.6 mg/mL extract load.