

For Research Use Only

Mouse KIM-1/HAVCR1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83064-1-PBS

Basic Information

Catalog Number: 83064-1-PBS	GenBank Accession Number: NM_134248	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 171283	CloneNo.: 230336A5
Source: Rabbit	UNIPROT ID: Q5QNS5	
Isotype: IgG	Full Name: hepatitis A virus cellular receptor 1	
	Calculated MW: 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Product Information

83064-1-PBS targets KIM-1/HAVCR1 as part of a matched antibody pair:

MP00116-2: 83064-1-PBS capture and 83064-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00116-2, Mouse KIM-1/HAVCR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83064-1-PBS. Detection antibody: 83064-2-PBS. Standard: Eg0571. Range: 0.625-20 ng/mL.