

For Research Use Only

SLC7A10 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83072-4-PBS

Basic Information

Catalog Number: 83072-4-PBS	GenBank Accession Number: BC035627	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 56301	CloneNo.: 230400G1
Source: Rabbit	UNIPROT ID: Q9NS82	
Isotype: IgG	Full Name: solute carrier family 7, (neutral amino acid transporter, y+ system) member 10	
Immunogen Catalog Number: AG34628	Calculated MW: 57 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83072-4-PBS targets SLC7A10 as part of a matched antibody pair:

MP00053-3: 83072-4-PBS capture and 83072-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00053-3, SLC7A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83072-4-PBS. Detection antibody: 83072-2-PBS. Standard: Ag34628. Range: 0.625-20 ng/mL.