

For Research Use Only

TNFR1/CD120a Recombinant antibody, PBS Only (Capture)

Catalog Number: 83080-1-PBS



Basic Information

Catalog Number: 83080-1-PBS	GenBank Accession Number: BC010140	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 7132	CloneNo.: 230322D3
Source: Rabbit	UNIPROT ID: P19438	
Isotype: IgG	Full Name: tumor necrosis factor receptor superfamily, member 1A	
	Calculated MW: 455 aa, 50 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83080-1-PBS targets TNFR1/CD120a as part of a matched antibody pair:

MP00099-1: 83080-1-PBS capture and 83080-3-PBS detection (validated in Cytometric bead array)

MP00099-2: 83080-1-PBS capture and 83080-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00099-1, TNFR1/CD120a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83080-1-PBS. Detection antibody: 83080-3-PBS. Standard: Eg0531. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00099-2, TNFR1/CD120a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83080-1-PBS. Detection antibody: 83080-2-PBS. Standard: Eg0531. Range: 1.25-20 ng/mL