

For Research Use Only

GDNF Recombinant antibody, PBS Only



Catalog Number: 83087-2-PBS

Basic Information

Catalog Number: 83087-2-PBS	GenBank Accession Number: BC069369	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2668	CloneNo.: 230358G6
Source: Rabbit	UNIPROT ID: P39905	
Isotype: IgG	Full Name: glial cell derived neurotrophic factor	
Immunogen Catalog Number: AG24068	Calculated MW: 211 aa, 24 kDa	

Applications

Tested Applications:
ELISA, WB

Species Specificity:
mouse, human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 83087-2-RR (GDNF antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83087-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83087-2-RR against Human GDNF were performed. The affinity constant is below 1 pM.