

For Research Use Only

DDX19A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83106-5-PBS

Basic Information

Catalog Number: 83106-5-PBS	GenBank Accession Number: BC005162	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 55308	CloneNo.: 230232B8
Source: Rabbit	UNIPROT ID: Q9NUU7	
Isotype: IgG	Full Name: DEAD (Asp-Glu-Ala-As) box polypeptide 19A	
Immunogen Catalog Number: AG7362	Calculated MW: 54 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83106-5-PBS targets DDX19A as part of a matched antibody pair:

MP00140-3: 83106-5-PBS capture and 83106-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00140-3, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-5-PBS. Detection antibody: 83106-4-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL.