

For Research Use Only

CHD6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83122-4-PBS

Basic Information

Catalog Number: 83122-4-PBS	GenBank Accession Number: NM_032221	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 84181	CloneNo.: 230397B12
Source: Rabbit	UNIPROT ID: Q8TD26	
Isotype: IgG	Full Name: chromodomain helicase DNA binding protein 6	
	Calculated MW: 305 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83122-4-PBS targets CHD6 as part of a matched antibody pair:

MP00298-2: 83122-6-PBS capture and 83122-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00298-2, CHD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83122-6-PBS. Detection antibody: 83122-4-PBS. Standard: SY01008. Range: 0.156-5 ng/mL