

For Research Use Only

RNF31 Recombinant antibody, PBS Only



Catalog Number: **83137-4-PBS**

Basic Information

Catalog Number: 83137-4-PBS	GenBank Accession Number: BC017376	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 55072	CloneNo.: 230446E9
Source: Rabbit	UNIPROT ID: Q96EPO	
Isotype: IgG	Full Name: ring finger protein 31	
Immunogen Catalog Number: AG34505	Calculated MW: 1072 aa, 120 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83137-4-RR (RNF31 / HOIP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83137-4-PBS in a different storage buffer formulation.