

For Research Use Only

CTNND2 Recombinant antibody, PBS Only



Catalog Number: 83138-5-PBS

Basic Information

Catalog Number:

83138-5-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001332

GeneID (NCBI):

1501

UNIPROT ID:

Q9UQB3

Full Name:

catenin (cadherin-associated protein), delta 2 (neural plakophilin-related arm-repeat protein)

Calculated MW:

133 kDa

Purification Method:

Protein A purification

CloneNo.:

230396F4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat brain tissue were subjected to SDS PAGE followed by western blot with 83138-5-RR (CTNND2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83138-5-PBS in a different storage buffer formulation.