

For Research Use Only

TIMP1 Recombinant antibody, PBS Only

Catalog Number: 83151-6-PBS



Basic Information

Catalog Number: 83151-6-PBS	GenBank Accession Number: BC000866	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7076	CloneNo.: 230329D8
Source: Rabbit	UNIPROT ID: P01033	
Isotype: IgG	Full Name: TIMP metalloproteinase inhibitor 1	
	Calculated MW: 23 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data