

For Research Use Only

FMO1 Recombinant antibody, PBS Only



Catalog Number: 83156-3-PBS

Basic Information

Catalog Number: 83156-3-PBS	GenBank Accession Number: NM_001282692	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 2326	CloneNo.: 230455B11
Source: Rabbit	UNIPROT ID: Q01740	
Isotype: IgG	Full Name: flavin containing monooxygenase 1	
Immunogen Catalog Number: AG34267	Calculated MW: 60 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83156-3-RR (FMO1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83156-3-PBS in a different storage buffer formulation.