

For Research Use Only

Mouse Angiopoietin-like 3/ANGPTL3 Recombinant antibody, PBS Only



Catalog Number: 83165-3-PBS

Basic Information

Catalog Number:

83165-3-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_013913.4

GeneID (NCBI):

30924

UNIPROT ID:

Q9R182

Full Name:

angiopoietin-like 3

Calculated MW:

53KD

Purification Method:

Protein A purification

CloneNo.:

230498F11

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83165-3-RR (Angiopoietin-like 3/ANGPTL3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83165-3-PBS in a different storage buffer formulation.



A375 cells were subjected to SDS PAGE followed by western blot with 83165-3-RR (Angiopoietin-like 3/ANGPTL3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83165-3-PBS in a different storage buffer formulation.