

For Research Use Only

PRKD1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83174-4-PBS

Basic Information

Catalog Number: 83174-4-PBS	GenBank Accession Number: NM_002742	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 5587	CloneNo.: 230355H3
Source: Rabbit	UNIPROT ID: Q15139	
Isotype: IgG	Full Name: protein kinase D1	
	Calculated MW: 102 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83174-4-PBS targets PRKD1 as part of a matched antibody pair:

MP00230-1: 83174-4-PBS capture and 83174-1-PBS detection (validated in Cytometric bead array)

MP00230-2: 83174-4-PBS capture and 83174-2-PBS detection (validated in Cytometric bead array)

MP00230-3: 83174-3-PBS capture and 83174-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve ofMP00230-1, PRKD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83174-4-PBS. Detection antibody: 83174-1-PBS. Standard: SY00930. Range: 0.156-2.5 ng/mL



Cytometric bead array standard curve ofMP00230-2, PRKD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83174-4-PBS. Detection antibody: 83174-2-PBS. Standard: SY00930. Range: 0.156-10 ng/mL



Cytometric bead array standard curve ofMP00230-3, PRKD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83174-3-PBS. Detection antibody: 83174-4-PBS. Standard: SY00930. Range: 0.156-10 ng/mL