

# TRIM54 Recombinant antibody, PBS Only (Capture)

Catalog Number: **83175-1-PBS**

## Basic Information

<b>Catalog Number:</b> 83175-1-PBS	<b>GenBank Accession Number:</b> BC141807	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 57159	<b>CloneNo.:</b> 230391A4
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9BYV2	
<b>Isotype:</b> IgG	<b>Full Name:</b> tripartite motif-containing 54	
<b>Immunogen Catalog Number:</b> AG15229	<b>Calculated MW:</b> 358 aa, 40 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Product Information

83175-1-PBS targets TRIM54 as part of a matched antibody pair:

MP00236-1: 83175-1-PBS capture and 83175-5-PBS detection (validated in Cytometric bead array)

MP00236-3: 83175-1-PBS capture and 83175-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00236-1, TRIM54 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83175-1-PBS. Detection antibody: 83175-5-PBS. Standard: Ag15229. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00236-3, TRIM54 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83175-1-PBS. Detection antibody: 83175-2-PBS. Standard: Ag15229. Range: 0.625-20 ng/mL