

For Research Use Only

MKS1 Recombinant antibody

Catalog Number: 83181-2-RR



Basic Information

Catalog Number: 83181-2-RR	GenBank Accession Number: BC010061	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 54903	CloneNo.: 240036E12
Source: Rabbit	UNIPROT ID: Q9NXB0	
Isotype: IgG	Full Name: Meckel syndrome, type 1	
Immunogen Catalog Number: AG9177	Calculated MW: 559 aa, 65 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83181-2-RR against Human BBS13/MKS1 were performed. The affinity constant is 0.198 nM.