

For Research Use Only

SLC38A9 Recombinant antibody, PBS Only



Catalog Number: 83182-5-PBS

Basic Information

Catalog Number:

83182-5-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15607

GenBank Accession Number:

BC066891

GeneID (NCBI):

153129

UNIPROT ID:

Q8NBW4

Full Name:

solute carrier family 38, member 9

Calculated MW:

561 aa, 64 kDa

Purification Method:

Protein A purification

CloneNo.:

240039E7

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with 83182-5-RR (SLC38A9 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83182-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83182-5-RR against Human SLC38A9 were performed. The affinity constant is 0.715 nM.