

For Research Use Only

PARK2/Parkin Recombinant antibody



Catalog Number: 83224-2-RR

Basic Information

Catalog Number: 83224-2-RR	GenBank Accession Number: BC022014	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 5071	CloneNo.: 240077C12
Source: Rabbit	UNIPROT ID: O60260	
Isotype: IgG	Full Name: Parkinson disease (autosomal recessive, juvenile) 2, parkin	
Immunogen Catalog Number: AG5179	Calculated MW: 52 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83224-2-RR against Human PARK2/Parkin were performed. The affinity constant is 89.4 pM.