

For Research Use Only

# SOX8 Recombinant antibody

Catalog Number: 83238-2-RR



## Basic Information

<b>Catalog Number:</b> 83238-2-RR	<b>GenBank Accession Number:</b> BC031797	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 30812	<b>CloneNo.:</b> 240120H5
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P57073	
<b>Isotype:</b> IgG	<b>Full Name:</b> SRY (sex determining region Y)-box 8	
<b>Immunogen Catalog Number:</b> AG34395	<b>Calculated MW:</b> 446 aa, 47 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83238-2-RR against Human SOX8 were performed. The affinity constant is 6.54 nM.