

For Research Use Only

SLC16A11 Recombinant antibody, PBS Only



Catalog Number: 83239-5-PBS

Basic Information

Catalog Number: 83239-5-PBS	GenBank Accession Number: BC093860	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 162515	CloneNo.: 240081H5
Source: Rabbit	UNIPROT ID: Q8NCK7	
Isotype: IgG	Full Name: solute carrier family 16, member 11 (monocarboxylic acid transporter 11)	
Immunogen Catalog Number: AG24920	Calculated MW: 471 aa, 48 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83239-5-RR (SLC16A11 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83239-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83239-5-RR against Human SLC16A11 were performed. The affinity constant is below 1 pM.