

For Research Use Only

NRG1, isoform Alpha Recombinant antibody, PBS Only (Detector)

Catalog Number: 83251-1-PBS



Basic Information

Catalog Number: 83251-1-PBS	GenBank Accession Number: BC007675	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 3084	CloneNo.: 240067A4
Source: Rabbit	UNIPROT ID: Q02297	
Isotype: IgG	Full Name: neuregulin 1	
Immunogen Catalog Number: AG35114	Calculated MW: 70 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

83251-1-PBS targets NRG1, isoform Alpha as part of a matched antibody pair:

MP00232-2: 83251-3-PBS capture and 83251-1-PBS detection (validated in Cytometric bead array)

MP00232-6: 83251-4-PBS capture and 83251-1-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP00232-6, Human NRG1, isoform Alpha Recombinant Matched Antibody Pair - PBS only. 83251-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag35114. 83251-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



The mean NRG1, isoform Alpha concentration was determined to be 2,318.1 pg/mL in A549 cell extract based on a 2.5 mg/mL extract load.



Cytometric bead array standard curve of MP00232-2, NRG1, isoform Alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83251-3-PBS. Detection antibody: 83251-1-PBS. Standard: Ag35114. Range: 1.25-20 ng/mL