

For Research Use Only

# DUS2L Recombinant antibody, PBS Only (Capture)

Catalog Number: 83265-2-PBS



## Basic Information

Catalog Number:

83265-2-PBS

Size:

100ug, Concentration: 1mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23831

GenBank Accession Number:

BC006527

GeneID (NCBI):

54920

UNIPROT ID:

Q9NX74

Full Name:

dihydrouridine synthase 2-like, SMM1  
homolog (S. cerevisiae)

Purification Method:

Protein A purification

CloneNo.:

240192F8

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83265-2-RR against Human DUS2L were performed. The affinity constant is 2.66 nM.