

For Research Use Only

DUS2L Recombinant antibody

Catalog Number: 83265-2-RR



Basic Information

Catalog Number: 83265-2-RR	GenBank Accession Number: BC006527	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 54920	CloneNo.: 240192F8
Source: Rabbit	UNIPROT ID: Q9NX74	
Isotype: IgG	Full Name: dihydrouridine synthase 2-like, SMM1 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG23831		

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83265-2-RR against Human DUS2L were performed. The affinity constant is 2.66 nM.