

For Research Use Only

TFPI2 Recombinant antibody, PBS Only



Catalog Number: 83279-1-PBS

Basic Information

Catalog Number: 83279-1-PBS	GenBank Accession Number: NM_001271003.1	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 7980	CloneNo.: 230135E3
Source: Rabbit	UNIPROT ID: P48307	
Isotype: IgG	Full Name: tissue factor pathway inhibitor 2	
Immunogen Catalog Number: AG33467	Calculated MW: 27 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 83279-1-RR (TFPI2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83279-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83279-1-RR against Human TFPI2 were performed. The affinity constant is below 1 pM.