

For Research Use Only

GAPDHS Recombinant antibody, PBS Only



Catalog Number: 83290-3-PBS

Basic Information

Catalog Number: 83290-3-PBS	GenBank Accession Number: BC036373	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 26330	CloneNo.: 240277F1
Source: Rabbit	UNIPROT ID: O14556	
Isotype: IgG	Full Name: glyceraldehyde-3-phosphate dehydrogenase, spermatogenic	
Immunogen Catalog Number: AG4988	Calculated MW: 45 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83290-3-RR (GAPDH antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83290-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83290-3-RR against Human GAPDH were performed. The affinity constant is below 1 pM.