

For Research Use Only

# GINS1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83293-3-PBS

## Basic Information

<b>Catalog Number:</b> 83293-3-PBS	<b>GenBank Accession Number:</b> BC012542	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, Concentration: 1mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 9837	<b>CloneNo.:</b> 240234B2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q14691	
<b>Isotype:</b> IgG	<b>Full Name:</b> GINS complex subunit 1 (Psf1 homolog)	
<b>Immunogen Catalog Number:</b> AG8895	<b>Calculated MW:</b> 196 aa, 23 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Product Information

83293-3-PBS targets GINS1 as part of a matched antibody pair:

MP00317-1: 83293-3-PBS capture and 83293-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)      W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00317-1, GINS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83293-3-PBS. Detection antibody: 83293-4-PBS. Standard: Ag8895. Range: 3.13-100 ng/mL