

For Research Use Only

APIP Recombinant antibody, PBS Only



Catalog Number: 83320-5-PBS

Basic Information

Catalog Number: 83320-5-PBS	GenBank Accession Number: BC009077	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51074	CloneNo.: 240189E2
Source: Rabbit	UNIPROT ID: Q96GX9	
Isotype: IgG	Full Name: APAF1 interacting protein	
Immunogen Catalog Number: AG2729	Calculated MW: 242 aa, 27 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83320-5-RR (AP1P antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83320-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83320-5-RR against Human AP1P were performed. The affinity constant is 5.86 nM.