

For Research Use Only

SLC38A7 Recombinant antibody, PBS Only



Catalog Number: 83346-6-PBS

Basic Information

Catalog Number:

83346-6-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35081

GenBank Accession Number:

BC001961

GeneID (NCBI):

55238

UNIPROT ID:

Q9NVC3

Full Name:

solute carrier family 38, member 7

Calculated MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

240147C3

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Caco-2 cells were subjected to SDS PAGE followed by western blot with 83346-6-RR (SLC38A7 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83346-6-PBS in a different storage buffer formulation.



human liver tissue were subjected to SDS PAGE followed by western blot with 83346-6-RR (SLC38A7 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83346-6-PBS in a different storage buffer formulation.