

For Research Use Only

CACHD1 Recombinant antibody, PBS Only



Catalog Number: 83350-4-PBS

Basic Information

Catalog Number: 83350-4-PBS	GenBank Accession Number: BC039301	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 57685	CloneNo.: 240197A2
Source: Rabbit	UNIPROT ID: Q5VU97	
Isotype: IgG	Full Name: cache domain containing 1	
Immunogen Catalog Number: AG22138	Calculated MW: 1274 aa, 142 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83350-4-RR (CACHD1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83350-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83350-4-RR against Human CACHD1 were performed. The affinity constant is 6.83 nM.