

For Research Use Only

TCIRG1 Recombinant antibody, PBS Only



Catalog Number: 83351-6-PBS

Basic Information

Catalog Number:

83351-6-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3378

GenBank Accession Number:

BC018133

GeneID (NCBI):

10312

UNIPROT ID:

Q13488

Full Name:

T-cell, immune regulator 1, ATPase, H+ transporting, lysosomal V0 subunit A3

Calculated MW:

830 aa, 93 kDa

Purification Method:

Protein A purification

CloneNo.:

240166B2

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83351-6-RR (TCIRG1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83351-6-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83351-6-RR against Human TCIRG1 were performed. The affinity constant is 2.14 nM.