

For Research Use Only

SLC22A8 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83357-4-PBS

Basic Information

| | | |
|---|---|---|
| Catalog Number: 83357-4-PBS | GenBank Accession Number: bc022387 | Purification Method: Protein A purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 9376 | CloneNo.: 240138G12 |
| Source: Rabbit | UNIPROT ID: Q8TCC7 | |
| Isotype: IgG | Full Name: solute carrier family 22 (organic anion transporter), member 8 | |
| Immunogen Catalog Number: AG9527 | Calculated MW: 542 aa, 60 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83357-4-PBS targets SLC22A8 as part of a matched antibody pair:

MP00363-1: 83357-4-PBS capture and 83357-5-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00363-1, SLC22A8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83357-4-PBS. Detection antibody: 83357-5-PBS. Standard: Ag9527. Range: 0.313-40ng/mL.