

For Research Use Only

# Cytokeratin 14 Recombinant antibody, PBS Only



Catalog Number: 83379-1-PBS

## Basic Information

Catalog Number:

83379-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0188

GenBank Accession Number:

BC002690

GeneID (NCBI):

3861

UNIPROT ID:

P02533

Full Name:

keratin 14

Calculated MW:

472 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

1H22

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



rat skin tissue were subjected to SDS PAGE followed by western blot with 83379-1-RR (KRT14 antibody) at dilution of 1:1400 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83379-1-PBS in a different storage buffer formulation.