

For Research Use Only

# SLC25A17 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83404-3-PBS

## Basic Information

**Catalog Number:**

83404-3-PBS

**Size:**

100ug, Concentration: 1 mg/ml by  
Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG34310

**GenBank Accession Number:**

BC012998

**GeneID (NCBI):**

10478

**UNIPROT ID:**

O43808

**Full Name:**

solute carrier family 25  
(mitochondrial carrier; peroxisomal  
membrane protein, 34kDa), member  
17

**Calculated MW:**

307 aa, 35 kDa

**Purification Method:**

Protein A purification

**CloneNo.:**

240382F5

## Applications

**Tested Applications:**

Indirect ELISA, Cytometric bead array

**Species Specificity:**

Human

## Product Information

83404-3-PBS targets SLC25A17 as part of a matched antibody pair:

MP00409-2: 83404-3-PBS capture and 83404-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00409-2, SLC25A17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83404-3-PBS. Detection antibody: 83404-2-PBS. Standard: Ag34310. Range: 2.5-80 ng/mL.