

For Research Use Only

TBX20 Recombinant antibody, PBS Only



Catalog Number: 83414-5-PBS

Basic Information

Catalog Number: 83414-5-PBS	GenBank Accession Number: BC120946	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 57057	CloneNo.: 240385C12
Source: Rabbit	UNIPROT ID: Q9UMR3	
Isotype: IgG	Full Name: T-box 20	
Immunogen Catalog Number: AG18501	Calculated MW: 297 aa, 33 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83414-5-RR (TBX20 antibody) at dilution of 1:15000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83414-5-PBS in a different storage buffer formulation.