

For Research Use Only

THOC5 Recombinant antibody

Catalog Number: 83435-1-RR



Basic Information

Catalog Number: 83435-1-RR	GenBank Accession Number: BC003615	Purification Method: Protein A purification
Size: 100ul , Concentration: 800 ug/ml by Nanodrop;	GeneID (NCBI): 8563	CloneNo.: 240013F7
Source: Rabbit	UNIPROT ID: Q13769	Recommended Dilutions: IF/ICC 1:200-1:800
Isotype: IgG	Full Name: THO complex 5	
Immunogen Catalog Number: AG6664	Calculated MW: 79 kDa	

Applications

Tested Applications: IF/ICC, FC (Intra), ELISA	Positive Controls: IF/ICC : U2OS cells,
Species Specificity: human	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ U2OS cells were intracellularly stained with 0.25 ug Thoc5 Recombinant Antibody (83435-1-RR, Clone:240013F7) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using THOC5 antibody (83435-1-RR, Clone: 240013F7) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using THOC5 antibody (83435-1-RR, Clone: 240013F7) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).