

For Research Use Only

GRIN2A Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83465-4-PBS



Basic Information

Catalog Number: 83465-4-PBS	GenBank Accession Number: NM_000833	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2903	CloneNo.: 240445F12
Source: Rabbit	UNIPROT ID: Q12879	
Isotype: IgG	Full Name: glutamate receptor, ionotropic, N-methyl D-aspartate 2A	
Immunogen Catalog Number: AG29101	Calculated MW: 165 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83465-4-PBS targets GRIN2A as part of a matched antibody pair:

MP00461-1: 83465-3-PBS capture and 83465-4-PBS detection (validated in Cytometric bead array)

MP00461-2: 83465-4-PBS capture and 83465-2-PBS detection (validated in Cytometric bead array)

MP00461-3: 83465-4-PBS capture and 83465-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Cytometric bead array standard curve of MP00461-1, GRIN2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83465-3-PBS. Detection antibody: 83465-4-PBS. Standard: Ag29101. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00461-2, GRIN2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83465-4-PBS. Detection antibody: 83465-2-PBS. Standard: Ag29101. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00461-3, GRIN2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83465-4-PBS. Detection antibody: 83465-1-PBS. Standard: Ag29101. Range: 0.313-40 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83465-4-PBS against Human GRIN2A was performed. The affinity constant is 3.99 nM.