

For Research Use Only

GSTT2B Recombinant antibody

Catalog Number: 83469-4-RR



Basic Information

Catalog Number: 83469-4-RR	GenBank Accession Number: BC071700	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 653689	CloneNo.: 240493G4
Source: Rabbit	UNIPROT ID: POCG30	Recommended Dilutions: WB 1:5000-1:50000
Isotype: IgG	Full Name: glutathione S-transferase theta 2B (gene/pseudogene)	
Immunogen Catalog Number: AG11828	Calculated MW: 244 aa, 28 kDa	
	Observed MW: 28 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human, mouse, rat

Positive Controls:
WB : HeLa cells, mouse liver tissue, mouse lung tissue,
rat liver tissue

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83469-4-RR (GSTT2B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 83469-4-RR against Human GSTT2B were performed. The affinity constant is below 1 pM.