

For Research Use Only

Mouse PD-1/CD279 Recombinant antibody, PBS Only



Catalog Number: 83516-6-PBS

Basic Information

Catalog Number: 83516-6-PBS	GenBank Accession Number: NM_008798.2	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 18566	CloneNo.: 240390D6
Source: Rabbit	UNIPROT ID: Q02242	
Isotype: IgG	Full Name: programmed cell death 1	
	Calculated MW: 32 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
mouse

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83516-6-RR (PD-1/CD279 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83516-6-PBS in a different storage buffer formulation.