

For Research Use Only

CD86 Recombinant antibody, PBS Only



Catalog Number: 83523-4-PBS

Basic Information

Catalog Number: 83523-4-PBS	GenBank Accession Number: NM_175862	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 942	CloneNo.: 240429G12
Source: Rabbit	ENSEMBL Gene ID: ENSG00000114013	
Isotype: IgG	UNIPROT ID: P42081	
	Full Name: CD86 molecule	
	Calculated MW: 329 aa, 38 kDa	

Applications

Tested Applications:
ELISA, WB

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Raji cells were subjected to SDS PAGE followed by western blot with 83523-4-RR (CD86 antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83523-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83523-4-RR against Human CD86 were performed. The affinity constant is 0.766 nM.