

For Research Use Only

METTL9 Recombinant antibody, PBS Only



Catalog Number: 83531-5-PBS

Basic Information

Catalog Number: 83531-5-PBS	GenBank Accession Number: BC000195	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51108	CloneNo.: 240510A2
Source: Rabbit	UNIPROT ID: Q9H1A3	
Isotype: IgG	Full Name: methyltransferase like 9	
Immunogen Catalog Number: AG7191	Calculated MW: 37 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 83531-5-RR (METTL9 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83531-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83531-5-RR against Human METTL9 were performed. The affinity constant is below 1 pM.