

For Research Use Only

SYTL4 Recombinant antibody, PBS Only



Catalog Number: 83564-5-PBS

Basic Information

Catalog Number: 83564-5-PBS	GenBank Accession Number: BC014913	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 94121	CloneNo.: 240556E2
Source: Rabbit	UNIPROT ID: Q96C24	
Isotype: IgG	Full Name: synaptotagmin-like 4	
Immunogen Catalog Number: AG2772	Calculated MW: 671 aa, 76 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83564-5-RR (SYTL4 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83564-5-PBS in a different storage buffer formulation.



Biacore BLI kinetic assays of 83564-5-RR against Human SYTL4 were performed. The affinity constant is below 1 pM.