

For Research Use Only

TRIM15 Recombinant antibody, PBS Only



Catalog Number: 83582-5-PBS

Basic Information

Catalog Number: 83582-5-PBS	GenBank Accession Number: BC038585	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 89870	CloneNo.: 240565C7
Source: Rabbit	UNIPROT ID: Q9C019	
Isotype: IgG	Full Name: tripartite motif-containing 15	
Immunogen Catalog Number: AG4533	Calculated MW: 465 aa, 52 kDa	

Applications

Tested Applications:
IF/ICC, ELISA

Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using TRIM15 antibody (83582-5-RR, Clone: 240565C7) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1). This data was developed using the same antibody clone with 83582-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83582-5-RR against Human TRIM15 were performed. The affinity constant is 0.686 nM.