

For Research Use Only

LEFTY1 Recombinant antibody, PBS Only



Catalog Number: 83586-1-PBS

Basic Information

Catalog Number: 83586-1-PBS	GenBank Accession Number: BC027883	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10637	CloneNo.: 240346E11
Source: Rabbit	UNIPROT ID: O75610	
Isotype: IgG	Full Name: left-right determination factor 1	
Immunogen Catalog Number: AG4925	Calculated MW: 366 aa, 41 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83586-1-RR (LEFTY1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83586-1-PBS in a different storage buffer formulation.