

For Research Use Only

CRY1 Recombinant antibody, PBS Only



Catalog Number: 83602-5-PBS

Basic Information

Catalog Number: 83602-5-PBS	GenBank Accession Number: BC030519	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1407	CloneNo.: 240520B12
Source: Rabbit	UNIPROT ID: Q16526	
Isotype: IgG	Full Name: cryptochrome 1 (photolyase-like)	
Immunogen Catalog Number: AG4270	Calculated MW: 586 aa, 66 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83602-5-RR (CRY1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83602-5-PBS in a different storage buffer formulation.