

For Research Use Only

# STEAP2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83611-4-PBS



## Basic Information

**Catalog Number:**

83611-4-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_152999

**GeneID (NCBI):**

261729

**UNIPROT ID:**

Q8NFT2

**Full Name:**

six transmembrane epithelial antigen of the prostate 2

**Calculated MW:**

56 kDa

**Purification Method:**

Protein A purification

**CloneNo.:**

240611F4

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Product Information

83611-4-PBS targets STEAP2 as part of a matched antibody pair:

MP00600-1: 83611-4-PBS capture and 83611-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00600-1, STEAP2-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83611-4-PBS. Detection antibody: 83611-2-PBS. Standard: SY00725. Range: 0.039-10 ng/mL.