

For Research Use Only

AP1M1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83613-4-PBS

Basic Information

Catalog Number: 83613-4-PBS	GenBank Accession Number: BC017469	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 8907	CloneNo.: 240553C1
Source: Rabbit	UNIPROT ID: Q9BXS5	
Isotype: IgG	Full Name: adaptor-related protein complex 1, mu 1 subunit	
Immunogen Catalog Number: AG2757	Calculated MW: 423 aa, 49 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83613-4-PBS targets AP1M1 as part of a matched antibody pair.

MP00588-2: 83613-2-PBS capture and 83613-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00588-2, AP1M1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83613-2-PBS. Detection antibody: 83613-4-PBS. Standard: Ag2757. Range: 3.125-100 ng/mL.