

For Research Use Only

CUL3 Recombinant antibody, PBS Only



Catalog Number:83619-5-PBS

Basic Information

Catalog Number: 83619-5-PBS	GenBank Accession Number: BC039598	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 8452	CloneNo.: 240534C7
Source: Rabbit	UNIPROT ID: Q13618	
Isotype: IgG	Full Name: cullin 3	
Immunogen Catalog Number: AG1555	Calculated MW: 89 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83619-5-RR (CUL3 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83619-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83619-5-RR against Human CUL3 were performed. The affinity constant is below 1 pM.