

For Research Use Only

# CDK8 Recombinant antibody, PBS Only

Catalog Number: 83621-6-PBS



## Basic Information

Catalog Number:

83621-6-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17148

GenBank Accession Number:

BC069634

GeneID (NCBI):

1024

UNIPROT ID:

P49336

Full Name:

cyclin-dependent kinase 8

Calculated MW:

464 aa, 53 kDa

Observed MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

240598G9

## Applications

Tested Applications:

WB, IF/ICC, FC (Intra), ELISA

Species Specificity:

human, mouse

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)


W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 83621-6-RR (CDK8 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83621-6-PBS in a different storage buffer formulation.



1x10<sup>6</sup> Caco-2 cells were intracellularly stained with 0.25 ug CDK8 Recombinant antibody (83621-6-RR, Clone:240598G9) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011). This data was developed using the same antibody clone with 83621-6-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CDK8 antibody (83621-6-RR, Clone: 240598G9) at dilution of 1:800 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83621-6-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83621-6-RR against Human CDK8 were performed. The affinity constant is 12.5 nM.