

For Research Use Only

CDK2 Recombinant antibody, PBS Only



Catalog Number: 83635-4-PBS

Basic Information

Catalog Number: 83635-4-PBS	GenBank Accession Number: BC003065	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1017	CloneNo.: 240580A11
Source: Rabbit	UNIPROT ID: P24941	
Isotype: IgG	Full Name: cyclin-dependent kinase 2	
Immunogen Catalog Number: AG17133	Calculated MW: 298 aa, 34 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, rat


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83635-4-RR (CDK2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83635-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83635-4-RR against Human CDK2 were performed. The affinity constant is 8.32 nM.