

For Research Use Only

RNASET2 Recombinant antibody, PBS Only



Catalog Number:83650-3-PBS

Basic Information

Catalog Number: 83650-3-PBS	GenBank Accession Number: BC039713	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 8635	CloneNo.: 240581H5
Source: Rabbit	UNIPROT ID: O00584	
Isotype: IgG	Full Name: ribonuclease T2	
Immunogen Catalog Number: AG28353	Calculated MW: 256 aa, 29 kDa	

Applications

Tested Applications:
ELISA, IF/ICC
Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed BxPC-3 cells using RNA5ET2 antibody (83650-3-RR, Clone: 240581H5) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83650-3-PBS in a different storage buffer formulation.